

Superabsorbent Dressing

Neosorb Superabsorbent Dressings quickly absorb and bind large amounts of wound exudate into the heart of the superabsorbent polymers. The outer envelope is ultrasonically bonded (adhesive-free construction) increasing conformability and minimising risk of possible skin reactions.



Ultrasonically Bonded

Adhesive and latex free, minimising risk of possible skin reactions



Highly Absorbent

Absorbs and retains wound exudate, even under compression



Sideways Swelling

Fluid uptake wicks sideways. Minimising change in the thickness of the dressing to ensure it is suitable for use under compression bandaging



Dry Contact Surface

Free draining wound contact surface facilitates fast absorption of fluid, that is then trapped within the superabsorbent polymer core



Suitable for:

- Moderate to highly exuding wounds
- Exudate management under compression therapy

| REF | Size | Sterility | Pcs/Unit | Units/Ctn |
|---------------|-------------|-----------|----------|-----------|
| NWC040 | 10cm x 10cm | Sterile | 20 | 20 |
| NWC046 | 10cm x 20cm | Sterile | 20 | 20 |
| NWC041 | 15cm x 15cm | Sterile | 20 | 20 |
| NWC042 | 15cm x 20cm | Sterile | 20 | 20 |
| NWC043 | 20cm x 20cm | Sterile | 20 | 20 |
| NWC044 | 20cm x 30cm | Sterile | 20 | 10 |
| NWC045 | 20cm x 45cm | Sterile | 5 | 20 |
| NWC047 | 30cm x 45cm | Sterile | 5 | 10 |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

