

# Liquid Control Dye Tests

Interster Liquid Control Dye Penetration Tests are specially developed to check the integrity of the seals produced with different heat sealers. Quality dye tests, as defined in ISO 11607-2, can be safely carried out without damaging the seal. Other advantages include being able to prepare and seal the test packaging in advance without the risk of leakage or migration of blue or red liquid on valuable objects in the work area.



**No Pipettes Used**  
Safely perform tests without damaging the seal



**Extended Shelf Life**  
5 years shelf life



**Suitable for:**

- **Type I:** All types of laminate packaging made of a film/medical grade paper combination
- **Type II:** All types of laminate packaging made of a film/Tyvek® combination

| REF         | Type                       | Ink Colour | Sterility | Pcs/Box | Boxes/Ctn |
|-------------|----------------------------|------------|-----------|---------|-----------|
| 3FSZB830816 | Type I: Film/Medical Paper | Blue       | Sterile   | 40      | 20        |
| 3FSZB830818 | Type II: Film/Tyvek®       | Red        | Sterile   | 40      | 20        |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.