

Lymphlex® REDUCE Reusable 2 Layer Compression Therapy Kit (2 Sets)



Reusable



Compress
& Supports



Reduces
Oedema

A convenient combination of reusable components offering a cost effective and environmentally responsible system to deliver two layer compression therapy of 30-40mm Hg for the management of Venous leg ulcers (VLUs).

REDUCE OEDEMA

Can be worn for up to 7 days, the clinician can choose how frequently compression is reapplied, dependent upon wound condition and oedema reduction

REDUCE WASTE

Independently tested for laundering up to 50 times - confidently reapply, as required, without increasing resources used or waste produced

REDUCE COST

Removes the need for single use compression bandaging, saving health facilities or the patient from unnecessary financial burden



Suitable for:

Venous leg ulcer management

Chronic oedema

Treatment of phlebological or lymphological conditions requiring compression

| Reorder | Components | Pcs/Unit | Units/Ctn |
|---------|--|----------|-----------|
| LOR100 | 2x Lymphlex AIR Reusable Foam Padding Bandage 8cm x 2.5m (LOR021) | 1 | 4 |
| | 2x Lymphlex AIR Reusable Foam Padding Bandage 10cm x 2.5m (LOR022) | | |
| | 2x Lymphlex Short Stretch Compression Bandage 8cm x 5m (LOR003) | | |
| | 2x Lymphlex Short Stretch Compression Bandage 10cm x 5m (LOR004) | | |
| | 1x Lymtube Tubular Overlay 10cm x 20cm (LOR043) | | |
| | 1x CoWrap Lite Lightweight Cohesive Bandage 8cm x 20m (LOR055) | | |
| | 1x Lymfix Fixation Tape 2.5cm x 5m (LOR050) | | |
| | 1x Laundry Bag | | |