

Lymphlex® TLC

Two Layer Cohesive Short Stretch (Inelastic) Compression Therapy



Easy Application



Compress & Supports



Low Profile

Lymphlex TLC Short Stretch Cohesive Kit combines the skin friendly Lymphlex AIR+ padding layer and Lymphlex SU short stretch cohesive compression bandages to create a safe, and easy to apply, compression therapy kit for the short term management of Venous Leg Ulcers. It is the perfect compression therapy option where laundering is problematic, or outside of local infection control guidelines.

THERAPEUTIC

Designed to provide a sub-bandage pressure of 30-40mmHg, when applied at full stretch, delivering consistent therapeutic compression

LOW PROFILE

Application of the 8cm bandage to the foot and ankle decreases the need to pleat and fold, creating a smooth and flat profile around this key part of the anatomy

COMFORT

Integrated skin protection layer, ensures even distribution of pressure whilst minimising skin irritation and reactivity, providing patient comfort in all seasons



Suitable for:

Venous leg ulcer management

Treatment of phlebological or lymphological conditions requiring compression

Reorder	Components	Pcs/Unit	Units/Ctn
LOS100 (LARGE)	1x Lymphlex AIR+ Padding Bandage 10cm x 4m (LOS022) 1x Lymphlex SU Short Stretch (Inelastic) Cohesive Compression Bandage 8cm x 6m (LOS003) 1x Lymphlex SU Short Stretch (Inelastic) Cohesive Compression Bandage 10cm x 6m (LOS004)	1	20
LOS101 (REGULAR)	1x Lymphlex AIR+ Padding Bandage 8cm x 4.5m (LOS021) 1x Lymphlex SU Short Stretch (Inelastic) Cohesive Compression Bandage 8cm x 6m (LOS003) 1x Lymphlex SU Short Stretch (Inelastic) Cohesive Compression Bandage 10cm x 6m (LOS004)	1	20