

Bolsaplast® Self Sealable Sterilisation Pouches

bolsaplast



Strong and
Durable



Various
Sizes



Self
Sealable

Produced with medical grade paper on one side and a PET/PP blend film on the other, the Bolsaplast Self Sealable Sterilisation Pouches offers extremely low fibre shedding and are resistant to film tearing. The pouches are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.

NON-TOXIC INK

No damage to instruments

COLOUR CHANGE INDICATOR

Chemical indicators for Steam/EO/Formaldehyde

MEDICAL GRADE CONSTRUCTION

Offers excellent peelability and is resistant to film tearing



Suitable for:

Sterilising devices and equipment

Reorder	Size	Pcs/Unit	Units/Ctn
MBA1080 <<	7.5cm x 20cm	200	6
MBA1081	9cm x 20cm	200	6
MBA1082	9cm x 23cm	200	6
MBA1084 <<	10cm x 24cm	200	6
MBA1085	13cm x 25cm	200	6
MBA1086	13cm x 33cm	200	6
MBA1087 <<	14cm x 28cm	200	4
MBA1088	20cm x 33cm	200	4
MBA1089	25cm x 40cm	200	4
MBA1090	30cm x 39cm	400	1
MBA1091 <<	30cm x 45cm	400	1

<< Indent Item

