



Strong and Durable



Various Sizes



Colour Change Indicator

Produced with Tyvek on one side and a PET/PP blend film on the other, Bolsaplast Heat Sealable Tyvek Pouches offer extremely low fibre shedding and are resistant to film tearing. The pouches are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.

COLOUR CHANGE INDICATOR

For plasma sterilisation (H₂O₂) printed under seal

NON-TOXIC INK

No damage to instruments

MEDICAL GRADE CONSTRUCTION

Offers excellent peelability and is resistant to film tearing



Suitable for:

Sterilising devices and equipment

Reorder	Size	Pcs/Unit	Units/Ctn
MBT1160<<	7.5cm x 21cm	1000	1
MBT1161<<	10cm x 28cm	1000	1
MBT1162<<	15cm x 32cm	1000	1
MBT1163<<	20cm x 42cm	350	1
MBT1164<<	25cm x 50cm	350	1
MBT1165<<	32cm x 56cm	400	1
MBT1166<<	40cm x 58cm	350	1

<< Indent Item