



Anti-Static

Low
LintingTear
Resistant

Senturian Minor Procedure Drape with Pigeonhole is designed to give a sterile field whilst providing access to the surgical site.

TEAR RESISTANT


Strong PE and non-woven composite fabric

MICROBIAL BARRIER

Decrease the risk of cross contamination and spread of hospital infections

CONVENIENT PIGEONHOLE

6.5cm pigeonhole with adhesive border



Suitable for use in procedure packs for:

Operating theatres

Emergency departments

Minor procedures kits

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
MDA001	50cm x 60cm	Non-Sterile	500	1