

## Sterilisation Flat Reel

Produced with medical grade paper on one side and a PET/PP blend film on the other, Bolsaplast Sterilisation Flat Reels offer extremely low fibre shedding and are resistant to film tearing. The reels are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.



### Colour Change Indicator

Chemical indicators for Steam/EO/formaldehyde



### Non-Toxic Ink

No damage to instruments



### Medical Grade Construction

Offers excellent peelability and is resistant to film tearing



### Compliance

Conforms to EN868-2/3/4/5 and ISO11607-1



**Suitable for:** - Sterilising devices and equipment

REF	Size	Sterility	Pcs/Unit	Units/Ctn
<b>MRP1040</b>	5cm x 200m	Non-sterile	1	1
<b>MRP1041</b>	7.5cm x 200m	Non-sterile	1	1
<b>MRP1042</b>	10cm x 200m	Non-sterile	1	1
<b>MRP1043</b> <<	12.5cm x 200m	Non-sterile	1	1
<b>MRP1044</b>	15cm x 200m	Non-sterile	1	1
<b>MRP1045</b>	20cm x 200m	Non-sterile	1	1
<b>MRP1046</b>	25cm x 200m	Non-sterile	1	1
<b>MRP1047</b>	30cm x 200m	Non-sterile	1	1
<b>MRP1048</b>	40cm x 200m	Non-sterile	1	1

<< Indent Item

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

