



Strong and
Durable



Various
Sizes



Colour Change
Indicator

Produced with Tyvek on one side and a PET/PP blend film on the other, Bolsaplast Flat Rolls offer extremely low fibre shedding and are resistant to film tearing. The rolls are manufactured under strict quality control guidelines to ensure consistent package integrity and indicator change after each sterilisation cycle.

COLOUR CHANGE INDICATOR

For plasma sterilisation (H₂O₂) printed under seal

NON-TOXIC INK

No damage to instruments

MEDICAL GRADE CONSTRUCTION

Offer excellent peelability and are resistant to film tearing



Suitable for:

Sterilising devices and equipment

Reorder	Size	Pcs/Unit	Units/Ctn
MRT1171	7.5cm x 70m	1	4
MRT1172	10cm x 70m	1	4
MRT1174	15cm x 70m	1	4
MRT1175	20cm x 70m	1	2
MRT1176	25cm x 70m	1	2
MRT1177<<	30cm x 70m	1	2
MRT1178<<	40cm x 70m	1	1
MRT1179<<	50cm x 70m	1	1

<< Indent Item