



Comfortable  
for Patient



Fast  
Absorbency



Gentle  
on Skin

Mediam 210 Low Adherent Dressing consists of an absorbent cotton inner pad encased between two layers of high density polyethylene film.

## LOW-ADHERENT

Low adherence to the wound, therefore minimal trauma and pain upon removal

## HIGHLY ABSORBENT

Highly absorbent inner pad, wicks excess exudate away from wound bed

## CUT TO SIZE

Maintains integrity when cut to size

## Suitable for:

Light to moderate exuding wounds

Cuts and grazes

Minor burns

Reorder	Size	Sterility	Pcs/Unit	Units/Ctn
NAD001	5cm x 5cm	Sterile	50	10
NAD002	7.5cm x 7.5cm	Sterile	50	10
NAD003	10cm x 10cm	Sterile	50	8
NAD004	10cm x 20cm	Sterile	50	8

