

# Sentry Medical Non-Woven Split Drainage Swab



Low Linting



Easy To Use



Highly Absorbent

Sentry Medical Non-Woven Split Drainage Swabs are constructed from a polyester and rayon blend. They have four times the absorbency of conventional cotton gauze and are low linting.

Cut with a Y split for easy placement around drains, they require no cutting and are ready to use.

## HIGHLY ABSORBENT

Provide four times the absorbency of conventional cotton gauze

## EASY TO USE

Y split for easy placement without cutting



## Suitable for:

Drain sites

Reorder	Pcs	Size	Sterility	Pcs/Unit	Units/Ctn
NWD1321	2pcs	7.5cm x 7.5cm	Sterile	50	10
NWD1421	2pcs	10cm x 10cm	Sterile	50	10

