

L2 Polybreath Impervious Gowns

Owear Polybreath Impervious Gowns are constructed from an impervious yet breathable material enhancing the wearer's comfort. These Polybreath Gowns meet the AAMI Level 2 liquid barrier performance levels, ensuring a reliable level of fluid resistance suitable for low to moderate risk scenarios.




AAMI Level 2

Gowns provide low barrier against fluid exposure and contaminants



Chemotherapy Tested

Provides protection when handling and administering chemotherapy*



Sonically Bonded

Strong sealed seams



Suitable for:

- Hospital and medical
- Light exposure to blood or body fluids
- Low risk patient isolation scenarios
- Pathology and ICU
- Laboratory

REF	Size	Feature	Colour	Pcs/Unit	Units/Ctn
GN009P	Regular	AAMI Level 2	White	50	1

*Suitable for use in preparing or handling chemotherapy drugs in oncology departments. Chemotherapy drugs tested against is available on request.

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

