

Level 2 Visor Mask

Anti-Fog - Anti-Glare - Procedure Ear Loops - Blue



PROTECTION STANDARDS



BFE (%):	≥ 98
PFE (%):	≥ 98
ΔP (mmH ₂ O/cm ²):	<6.0
Fluid Resistance:	120mm Hg
Flame Spread:	Class 1

Variation in lot results exist.
All results must meet or exceed
ASTM F2100 Requirements.

Based on ASTM F2100 standards, this mask provides a moderate level of protection during medical procedures. The attached visor and fluid-resistant and filtration capabilities of the mask provide all-in-one protection for wound care, line changes, dental work and procedures with a moderate risk of fluid splash when handling trauma or caring for ER or ICU patients.

SOFT MASK AND EAR LOOP MATERIAL

Carefully engineered and clinically preferred for extended periods of use without discomfort

ANTI-FOG VISOR

Designed with wrap-around wings and anti-fog properties, providing all-in-one splash protection with an unobstructed view

WIDE NOSEPIECE WITH ANTI-FOG FOAM

Constructed with high-quality aluminum, anti-fog foam and a black anti-glare strip, providing enhanced visibility and comfort

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin



Suitable for:

Hospital and medical

Emergency rooms

Areas with moderate risk of fluid splash

Reorder	Size	Feature	Pcs/Unit	Units/Ctn
PG4-1592	One Size	Level 2	25	4



Reorder	Description	Colour	Anti-fog Feature	Fluid Resistance ¹	BFE ²	Pcs/Box	Boxes/Case	Pcs/Case
LEVEL 2 BARRIER PROTECTION - FLUID RESISTANCE 120MM HG								
PG4-1592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Procedure Ear Loops	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PG4-2592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Surgical Ties	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PM4-305	ASTM Level 2 Mask / Surgical Ties	Blue	-	120mm Hg	≥ 98%	50	6	300
PM4-306	ASTM Level 2 Mask / Procedure Ear Loops	Blue	-	120mm Hg	≥ 98%	50	10	500
PM4-311	ASTM Level 2 Anti-Fog Strip Mask / Surgical Ties	Blue	Film	120mm Hg	≥ 98%	50	6	300
PM4-312	ASTM Level 2 Anti-Fog Strip Mask / Procedure Ear Loops	Blue	Film	120mm Hg	≥ 98%	50	10	500
LEVEL 3 BARRIER PROTECTION - FLUID RESISTANCE 160 MM HG								
PG4-1073	ASTM Level 3 Anti-fog Foam Mask / Procedure Ear Loops	Indigo	Foam	160mm Hg	≥ 98%	50	10	500
PG4-1273	ASTM Level 3 Mask / Procedure Ear Loops	Blue	-	160mm Hg	≥ 98%	50	10	500
PG4-1333	ASTM Level 3 Soft Mask / Procedure Ear Loops / Irritant-Free	White	-	160mm Hg	≥ 98%	50	10	500
PG4-1563	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Procedure Ear Loops	Blue	Foam	160mm Hg	≥ 98%	25	4	100
PG4-2333	ASTM Level 3 Soft Mask / Surgical Ties / Irritant-Free	White	-	160mm Hg	≥ 98%	50	6	300
PG4-3033	ASTM Level 3 Irritant-Free Duckbill Mask / Surgical Ties	White	-	160mm Hg	≥ 98%	50	6	300
PG4-4053	ASTM Level 3 Anti-Fog Tape Mask / Surgical Ties	Green	Adhesive Tape	160mm Hg	≥ 98%	50	6	300
PG4-4073	ASTM Level 3 Anti-fog Foam Mask / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	50	6	300
PG4-5073	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	25	4	100
PG4-6003	ASTM Level 3 High-Filtration Mask / Surgical Ties	Grey	-	160mm Hg	≥ 98%	50	6	300
PAEDIATRIC LINE								
PG4-1680	Paediatric Mask / Dinosaur Pattern	Dinosaurs	-	-	-	50	10	500

1. Resistance to Penetration by Synthetic Blood. ASTM 1862

2. Bacterial Filtration Efficiency. ASTM F2101