

Paediatric Mask Ear Loops - Dinosaur Pattern



Easy To Use



Suitable For Children



Breathable

A go-to mask for respiratory etiquette for symptomatic children who may be coughing or sneezing. This mask is specifically designed in a smaller size to ensure a proper fit on kids and the dinosaur pattern gives the mask a more friendly look to put young patients at ease.

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the skin

WIDE NOSEPIECE

Constructed with high-quality aluminum to keep the mask in place across the nose and cheekbones

SOFT MASK AND EAR LOOP MATERIAL

Carefully engineered for extended periods of use without discomfort



Suitable for:

Hospital and medical

Flu-like symptoms

Visiting hospital / medical environment

Reorder	Size	Ages	Pcs/Unit	Units/Ctn
PG4-1680	One Size	Children	50	10



Reorder	Description	Colour	Anti-fog Feature	Fluid Resistance ¹	BFE ²	Pcs/Box	Boxes/Case	Pcs/Case
LEVEL 2 BARRIER PROTECTION - FLUID RESISTANCE 120MM HG								
PG4-1592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Procedure Ear Loops	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PG4-2592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Surgical Ties	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PM4-305	ASTM Level 2 Mask / Surgical Ties	Blue	-	120mm Hg	≥ 98%	50	6	300
PM4-306	ASTM Level 2 Mask / Procedure Ear Loops	Blue	-	120mm Hg	≥ 98%	50	10	500
PM4-311	ASTM Level 2 Anti-Fog Strip Mask / Surgical Ties	Blue	Film	120mm Hg	≥ 98%	50	6	300
PM4-312	ASTM Level 2 Anti-Fog Strip Mask / Procedure Ear Loops	Blue	Film	120mm Hg	≥ 98%	50	10	500
LEVEL 3 BARRIER PROTECTION - FLUID RESISTANCE 160 MM HG								
PG4-1073	ASTM Level 3 Anti-fog Foam Mask / Procedure Ear Loops	Indigo	Foam	160mm Hg	≥ 98%	50	10	500
PG4-1273	ASTM Level 3 Mask / Procedure Ear Loops	Blue	-	160mm Hg	≥ 98%	50	10	500
PG4-1333	ASTM Level 3 Soft Mask / Procedure Ear Loops / Irritant-Free	White	-	160mm Hg	≥ 98%	50	10	500
PG4-1563	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Procedure Ear Loops	Blue	Foam	160mm Hg	≥ 98%	25	4	100
PG4-2333	ASTM Level 3 Soft Mask / Surgical Ties / Irritant-Free	White	-	160mm Hg	≥ 98%	50	6	300
PG4-3033	ASTM Level 3 Irritant-Free Duckbill Mask / Surgical Ties	White	-	160mm Hg	≥ 98%	50	6	300
PG4-4053	ASTM Level 3 Anti-Fog Tape Mask / Surgical Ties	Green	Adhesive Tape	160mm Hg	≥ 98%	50	6	300
PG4-4073	ASTM Level 3 Anti-fog Foam Mask / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	50	6	300
PG4-5073	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	25	4	100
PG4-6003	ASTM Level 3 High-Filtration Mask / Surgical Ties	Grey	-	160mm Hg	≥ 98%	50	6	300
PAEDIATRIC LINE								
PG4-1680	Paediatric Mask / Dinosaur Pattern	Dinosaurs	-	-	-	50	10	500

1. Resistance to Penetration by Synthetic Blood. ASTM 1862

2. Bacterial Filtration Efficiency. ASTM F2101