

# Level 2 Visor Mask

## Anti-Fog - Anti-Glare - Surgical Ties - Blue



### PROTECTION STANDARDS



BFE (%):  $\geq 98$   
 PFE (%):  $\geq 98$   
 $\Delta P$  (mmH<sub>2</sub>O/cm<sup>2</sup>):  $<6.0$   
 Fluid Resistance: 120mm Hg  
 Flame Spread: Class 1

Variation in lot results exist.  
 All results must meet or exceed  
 ASTM F2100 Requirements.

This mask provides protection during surgical procedures with moderate risk of fluid splash. It creates an important barrier for a patient against contaminants generated by clinicians and protects clinicians from fluids or particles from the patient. The attached visor provides additional anti-glare and anti-fog eye protection.

### LONG SURGICAL TIES

Long ties allow for easy donning and doffing, with flexibility to adjust the tightness and fit of the mask over a surgical or bouffant cap

### ANTI-FOG VISOR

Designed with wrap-around wings and anti-fog properties, providing all-in-one splash protection with an unobstructed view

### WIDE NOSEPIECE WITH ANTI-FOG FOAM

Constructed with high-quality aluminum, anti-fog foam and a black anti-glare strip, providing enhanced visibility and comfort

### ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin



### Suitable for:

Hospital and medical

Emergency rooms

Areas with moderate risk of fluid splash

Reorder	Size	Feature	Pcs/Unit	Units/Ctn
PG4-2592	One Size	Level 2	25	4



Reorder	Description	Colour	Anti-fog Feature	Fluid Resistance <sup>1</sup>	BFE <sup>2</sup>	Pcs/Box	Boxes/Case	Pcs/Case
<b>LEVEL 2 BARRIER PROTECTION - FLUID RESISTANCE 120MM HG</b>								
PG4-1592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Procedure Ear Loops	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PG4-2592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Surgical Ties	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PM4-305	ASTM Level 2 Mask / Surgical Ties	Blue	-	120mm Hg	≥ 98%	50	6	300
PM4-306	ASTM Level 2 Mask / Procedure Ear Loops	Blue	-	120mm Hg	≥ 98%	50	10	500
PM4-311	ASTM Level 2 Anti-Fog Strip Mask / Surgical Ties	Blue	Film	120mm Hg	≥ 98%	50	6	300
PM4-312	ASTM Level 2 Anti-Fog Strip Mask / Procedure Ear Loops	Blue	Film	120mm Hg	≥ 98%	50	10	500
<b>LEVEL 3 BARRIER PROTECTION - FLUID RESISTANCE 160 MM HG</b>								
PG4-1073	ASTM Level 3 Anti-fog Foam Mask / Procedure Ear Loops	Indigo	Foam	160mm Hg	≥ 98%	50	10	500
PG4-1273	ASTM Level 3 Mask / Procedure Ear Loops	Blue	-	160mm Hg	≥ 98%	50	10	500
PG4-1333	ASTM Level 3 Soft Mask / Procedure Ear Loops / Irritant-Free	White	-	160mm Hg	≥ 98%	50	10	500
PG4-1563	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Procedure Ear Loops	Blue	Foam	160mm Hg	≥ 98%	25	4	100
PG4-2333	ASTM Level 3 Soft Mask / Surgical Ties / Irritant-Free	White	-	160mm Hg	≥ 98%	50	6	300
PG4-3033	ASTM Level 3 Irritant-Free Duckbill Mask / Surgical Ties	White	-	160mm Hg	≥ 98%	50	6	300
PG4-4053	ASTM Level 3 Anti-Fog Tape Mask / Surgical Ties	Green	Adhesive Tape	160mm Hg	≥ 98%	50	6	300
PG4-4073	ASTM Level 3 Anti-fog Foam Mask / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	50	6	300
PG4-5073	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	25	4	100
PG4-6003	ASTM Level 3 High-Filtration Mask / Surgical Ties	Grey	-	160mm Hg	≥ 98%	50	6	300
<b>PAEDIATRIC LINE</b>								
PG4-1680	Paediatric Mask / Dinosaur Pattern	Dinosaurs	-	-	-	50	10	500

1. Resistance to Penetration by Synthetic Blood. ASTM 1862

2. Bacterial Filtration Efficiency. ASTM F2101