

Level 2 Anti-Fog Strip Mask Surgical Ties - Blue



PROTECTION STANDARDS



BFE (%):	≥ 98
PFE (%):	≥ 98
Δ P (mmH ₂ O/cm ²):	<6.0
Fluid Resistance:	120mm Hg
Flame Spread:	Class 1

Variation in lot results exist.
All results must meet or exceed
ASTM F2100 Requirements.

This mask provides a moderate level of protection during surgical procedures. It creates an important barrier for a patient against contaminants generated by clinicians and protects clinicians from fluids or particles from the patient. Gentle anti-fog film clings to the skin to trap humidity, ensuring a clear line of sight for clinicians wearing glasses or goggles during a procedure.

LONG SURGICAL TIES

Long ties allow for easy donning and doffing, with flexibility to adjust the tightness and fit of the mask over a surgical or bouffant cap

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin

WIDE NOSEPIECE WITH ANTI-FOG FILM

Constructed with high-quality aluminum and a gentle dual plastic anti-fog film designed to cling to the skin and deflect exhaled humidity away from the wearer's glasses



Suitable for:

Hospital and medical

Emergency rooms

Areas with moderate risk of fluid splash

Reorder	Size	Feature	Pcs/Unit	Units/Ctn
PM4-311	One Size	Level 2	50	6



Reorder	Description	Colour	Anti-fog Feature	Fluid Resistance ¹	BFE ²	Pcs/Box	Boxes/Case	Pcs/Case
LEVEL 2 BARRIER PROTECTION - FLUID RESISTANCE 120MM HG								
PG4-1592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Procedure Ear Loops	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PG4-2592	ASTM Level 2 Visor Mask / Anti-Fog / Anti-Glare / Surgical Ties	Blue	Foam	120mm Hg	≥ 98%	25	4	100
PM4-305	ASTM Level 2 Mask / Surgical Ties	Blue	-	120mm Hg	≥ 98%	50	6	300
PM4-306	ASTM Level 2 Mask / Procedure Ear Loops	Blue	-	120mm Hg	≥ 98%	50	10	500
PM4-311	ASTM Level 2 Anti-Fog Strip Mask / Surgical Ties	Blue	Film	120mm Hg	≥ 98%	50	6	300
PM4-312	ASTM Level 2 Anti-Fog Strip Mask / Procedure Ear Loops	Blue	Film	120mm Hg	≥ 98%	50	10	500
LEVEL 3 BARRIER PROTECTION - FLUID RESISTANCE 160 MM HG								
PG4-1073	ASTM Level 3 Anti-fog Foam Mask / Procedure Ear Loops	Indigo	Foam	160mm Hg	≥ 98%	50	10	500
PG4-1273	ASTM Level 3 Mask / Procedure Ear Loops	Blue	-	160mm Hg	≥ 98%	50	10	500
PG4-1333	ASTM Level 3 Soft Mask / Procedure Ear Loops / Irritant-Free	White	-	160mm Hg	≥ 98%	50	10	500
PG4-1563	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Procedure Ear Loops	Blue	Foam	160mm Hg	≥ 98%	25	4	100
PG4-2333	ASTM Level 3 Soft Mask / Surgical Ties / Irritant-Free	White	-	160mm Hg	≥ 98%	50	6	300
PG4-3033	ASTM Level 3 Irritant-Free Duckbill Mask / Surgical Ties	White	-	160mm Hg	≥ 98%	50	6	300
PG4-4053	ASTM Level 3 Anti-Fog Tape Mask / Surgical Ties	Green	Adhesive Tape	160mm Hg	≥ 98%	50	6	300
PG4-4073	ASTM Level 3 Anti-fog Foam Mask / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	50	6	300
PG4-5073	ASTM Level 3 Visor Mask / Anti-fog / Anti-glare / Surgical Ties	Indigo	Foam	160mm Hg	≥ 98%	25	4	100
PG4-6003	ASTM Level 3 High-Filtration Mask / Surgical Ties	Grey	-	160mm Hg	≥ 98%	50	6	300
PAEDIATRIC LINE								
PG4-1680	Paediatric Mask / Dinosaur Pattern	Dinosaurs	-	-	-	50	10	500

1. Resistance to Penetration by Synthetic Blood. ASTM 1862

2. Bacterial Filtration Efficiency. ASTM F2101