

Sentry Medical Resuscitation Face Shield and Mask



One Way
Valve



Ideal For
First Aid



Eliminates Cross
Contamination

Assists in the application of mouth to mouth resuscitation in emergency and first aid situations. These shields incorporate a one way valve to allow air to be transferred into the patient without cross contamination.

ONE WAY AIR VALVE

For patient and user safety

COMPACT AND LIGHTWEIGHT

Ideal for first aid and emergency kits

Suitable for:

First aid

Emergency



Reorder	Description	Size	Sterility	Pcs/Unit	Units/Ctn
RES001	Resuscitation Face Shield	One Size	Non-Sterile	50	12
10120245	S+M Resuscitation Mask with Oxygen Port	One Size	Non-Sterile	50	1

