



Sage[®] Anaesthetic Packs

Features

Sterility guaranteed

Quality components

No CSSD requirements

Complete one step pack

Small pack size – minimal storage

Puncture resistant tray with no pinholes

Convenient and easy to use, aseptic design

Standardisation of procedure requirements

Environmentally safe - incineration process converts procedure packs to harmless gasses of oxygen and carbon dioxide



PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
11002006	Sage [®] Anaesthetic Pack #1 Sterile	200	20	200
11002014	Sage [®] Anaesthetic #2 Sterile	75	N/A	75
11002022	Sage [®] Anaesthetic Pack #6 Sterile	100	N/A	100



Sage[®] Anaesthetic Pack

Anaesthetic Pack #6 - Tear Pack

PRODUCT CODE: 11002006

1 x Multipurpose Tray 220 x 120 x 20mm

2 x Green Gauze Swabs 10 x 10cm 8-Ply

1 x Towel 2-Ply Paper 42 x 40cm



Anaesthetic Dressing Pack #14 - Peel Pouch

PRODUCT CODE: 11002014

1 x Multipurpose Tray 220 x 120 x 20mm

5 x Green Gauze Swabs 10 x 10cm 8-Ply

1 x Forceps Green

1 x Drawing Up Needle 18g x 1.5" (38mm)

1 x Hypo Needle 19g x 1.5" (38mm)

1 x Hypo Needle 25g x 1" (25mm)

1 x Syringe 2/3ml Luer Slip

1 x Syringe 5ml Luer Slip

1 x Fenestrated Plastic Drape 70 x 100cm with Double Sided Tape Fenestration (13 x 13cm)

1 x Towel Scrim Reinforced 33 x 33cm

1 x Field White Plastic 75 x 75cm



Anaesthetic Dressing Pack #22 - Tear Pack

PRODUCT CODE: 11002022

1 x Multipurpose Tray 220 x 120 x 20mm

6 x Gauze Swabs Green 10 x 10cm

2 x Towel 2-Ply Paper 40 x 40cm

1 x Field White Plastic 50 x 50cm

