



Sage[®] Catheter Packs

Features

Complete one step tear pack

Puncture resistant tray

No CSSD requirements

Environmentally safe

Quality components

Sterility guaranteed

Increased efficiency

Easy to use

Cost effective



PRODUCT CODE	PRODUCT DESCRIPTION		SKU (MOQ)	SKU/ CTN
11003001	Sage [®] Catheter Pack #1	Sterile	30	30
11003002	Sage [®] Catheter Pack #2	Sterile	60	60
11003006	Sage [®] Catheter Pack #6	Sterile	60	60
11003008	Sage [®] Catheter Pack #8	Sterile	60	60
11003009	Sage [®] Catheter Pack #9	Sterile	55	55

Sage[®] Catheter Packs

Catheter Pack #1

PRODUCT CODE: 11003001

-
- 2 x Kidney Dish 375ml Green
-
- 2 x Gallipots 150ml
-
- 20 x Cotton Balls
-
- 2 x Gauze Swabs 10 x 10cm 8-Ply
-
- 4 x Forceps Green
-
- 2 x Towel 2-Ply Paper 42 x 40cm
-
- 1 x Field White Plastic 75 x 75cm
-



Catheter Pack #2

PRODUCT CODE: 11003002

-
- 2 x Long Tray 230 x 90 x 35mm
-
- 7 x Cotton Balls
-
- 2 x Forceps Green
-
- 2 x Examination Gloves Medium Nitrile Powder/Free
-
- 1 x Fenestrated Drape 4-Ply Polylined 51 x 68cm
(Fenestration 10 x 10cm)
-
- 1 x Towel 2-Ply Paper 42 x 40cm
-
- 1 x Field White Plastic 60 x 60cm
-



Catheter Pack #6

PRODUCT CODE: 11003006

-
- 1 x Long Tray 230 x 90 x 35mm
-
- 6 x Gauze Swabs 10 x 10cm 8-Ply
-
- 1 x Lube Jel 3gm
-
- 1 x Water Ampoule for Injections 10ml
-
- 1 x Saline Ampoule 30ml
-
- 1 x Syringe 10ml Eccentric Luer Slip
-
- 2 x Nitrile Gloves Large Powder Free
-
- 1 x Fenestrated Blue Drape Cellulose/Polyester 60 x 60cm
(Fenestration 10 x 10cm)
-
- 1 x Towel 2-Ply Paper 42 x 40cm
-
- 1 x Field White Plastic 60 x 60cm
-



Sage[®] Catheter Packs

Catheter Pack #8

PRODUCT CODE: 11003008

- 1 x Long Tray 230 x 90 x 35mm
- 5 x Cotton Balls
- 6 x Gauze Swabs 7.5 x 7.5cm 8-Ply
- 3 x Forceps Green
- 1 x Lube Jelly 2.7gm
- 1 x Water Ampoule for Injections 10ml
- 1 x Saline Ampoule 30ml
- 2 x Syringe 10ml Eccentric Luer Slip
- 2 x Examination Gloves Large Nitrile Powder-Free
- 1 x Fenestrated Blue Drape Cellulose/Polyester 60 x 60cm (Fenestration 10 x 10cm)
- 1 x Towel 2-Ply Paper 42 x 40cm
- 1 x Field White Plastic 60 x 60cm



Catheter Pack #9

PRODUCT CODE: 11003009

- 2 x Long Tray 230 x 90 x 35mm
- 5 x Gauze Swabs 7.5 x 7.5cm 8-Ply
- 2 x Forceps Green
- 1 x Forceps Yellow
- 1 x Fenestrated Drape 4-Ply Polylined 51 x 67cm (Fenestration 10 x 10cm)
- 1 x Towel 4-Ply Paper 42 x 40cm
- 1 x White Plastic Field 60 x 60cm

