



Sage[®] Lumbar Packs – Infant

Features

Peel pouch guarantees sterility

Variety of gauges across all our needles

Choice of Quinke, Hypo or Sprotte needles

Convenient kit form has all requirements and reduces preparation time



PRODUCT CODE	PRODUCT DESCRIPTION	SKU (MOQ)	SKU/ CTN
11015013	Sage [®] Lumbar Pack Infant #2 Sterile	10	10
11015014	Sage [®] Lumbar Pack infant #3 Sterile	10	10
11015016	Sage [®] Lumbar Pack Infant #5 Sterile	10	10
11015017	Sage [®] Lumbar Pack Infant #6 Sterile	10	10



Sage[®] Lumbar Packs – Infant

Lumbar Pack #2

PRODUCT CODE: 11015013

1 x Large Plain Tray 250 x 200 x 25mm
1 x Small (Non Sterile) Plain Tray 180 x 120 x 15mm
2 x Gallipots 60ml
5 x Gauze Swabs 5 x 5cm 8-Ply
1 x Spinal Needle 22g x 1.5" (38mm)
1 x Syringe 2 ml Luer Slip CN
1 x Hypo Needle 25g x 5/8" (16mm)
2 x Forceps Green
1 x Stopcock 3-Way White Tap Luer Lock
1 x Specitube 5ml Yellow Cap #1
1 x Specitube 10ml Yellow Cap #2
1 x Specitube 10ml Yellow Cap #3
1 x Strip Bandage 19 x 72mm
1 x Spinal Manometer
0.5 x Petri-Dish
10 x Cotton Balls 0.6gm
1 x Towel Paper 4-Ply 42 x 40cm
1 x Fenestrated Green Crepe Drape 60 x 30cm
(Fenestration 4 x 8cm)
1 x Field Green Crepe 60 x 60cm

Lumbar Pack #5

PRODUCT CODE: 11015016

1 x Large Plain Tray 250 x 200 x 25mm
1 x Small (Non Sterile) Plain Tray 180 x 120 x 15mm
2 x Gallipots 60ml
5 x Non-Woven Swabs 5 x 5 cm 4 Ply
1 x Sprotte Cannula 25G x 1 3/8" (35mm) (Latex Free)
2 x Forceps Green
1 x Syringe 2ml Luer Slip CN
1 x Hypo Needle 25g x 5/8" (16mm)
3 x Alcohol Swabs
1 x Stopcock 3-Way White Tap L/L
1 x Specitube 5ml Yellow Cap #1
1 x Specitube 10ml Yellow Cap #2
1 x Specitube 10ml Yellow Cap #3
1 x Specitube 5ml Yellow Cap #4
1 x Strip Bandage 19 x 72mm
1 x Spinal Manometer
0.5 x Petri-Dish
10 x Cotton Balls 0.6gm
1 x Towel 4-Ply Paper 42 x 40cm
1 x Plastic Drape 75 x 100cm D/S Tape (Fen 12.5 x 12.5cm)
1 x Field Green Crepe 60 x 60cm

Lumbar Pack #3

PRODUCT CODE: 11015014

1 x Large Plain Tray 250 x 200 x 25mm
1 x Small (Non Sterile) Plain Tray 180 x 120 x 15mm
2 x Gallipots 60ml
5 x Gauze Swabs 5 x 5 cm 8-Ply
1 x Sprotte cannula 25g x 1 3/8" (35mm)
2 x Forceps Green
1 x Syringe 2 ml Luer Slip CN
1 x Hypo Needle 25g x 5/8" (16mm)
1 x Stopcock 3-Way White Tap L/L
1 x Specitube 5ml Yellow Cap #1
1 x Specitube 10ml Yellow Cap #2
1 x Specitube 10ml Yellow Cap #3
1 x Strip Bandage 19 x 72mm
1 x Spinal Manometer
0.5 x Petri-Dish
10 x Cotton Balls 0.6gm
1 x Towel Paper 4-Ply 42 x 40cm
1 x A&B - Field Green Crepe 60 x 60cm

Lumbar Pack #6

PRODUCT CODE: 11015017

1 x Large Plain Tray 250 x 200 x 25mm
1 x Small (Non Sterile) Plain Tray 180 x 120 x 15mm
2 x Gallipot 60ml
5 x Non-Woven Swabs 5 x 5 cm 4-Ply
2 x Forceps Green
1 x Syringe 2ml L/Slip CN
1 x Hypo Needle 25g x 5/8" (16mm)
3 x Alcohol Swab
1 x Stopcock 3-Way White Tap L/L
1 x Specitube 5ml Yellow Cap #1
1 x Specitube 10ml Yellow Cap #2
1 x Specitube 10ml Yellow Cap #3
1 x Specitube 5ml Yellow Cap #4
1 x Strip Bandage 19 x 72mm
1 x Spinal Manometre
0.5 x Petri Dish
10 x Cotton Ball 0.6gm
1 x Towel Paper 4-Ply 42 x 40cm
1 x Plastic Drape 75 x 100cm D/S Tape (Fen 12.5 x 12.5cm)
1 x Field Green Crepe 60 x 60cm