



Sage[®] Suture Tray Packs

Features

No CSSD requirements

Peel pouch guarantees sterility

Small pack size – minimal storage

Puncture resistant tray with no pinholes

Convenient and easy to use, aseptic design

Environmentally safe - incineration process converts procedure packs to harmless gasses of oxygen and carbon dioxide

Standardisation of procedure requirements reduces preparation time



PRODUCT CODE	PRODUCT DESCRIPTION	SKU (MOQ)	SKU/ Inner	SKU/ CTN
11011049	Sage [®] Suture Pack #3 Sterile	10	N/A	10
11011052	Sage [®] Suture Pack #6 Sterile	10	N/A	10
11011057	Sage [®] Suture Pack Removal Sterile	120	N/A	120
11011062	Sage [®] Suture Pack #7 Sterile	10	N/A	10
11011131	Sage [®] Suture Pack #12 Sterile	10	N/A	10
11011133	Sage [®] Suture Pack Removal #1 Sterile	35	35	350
11011157	Sage [®] Suture Pack #13 Sterile	10	N/A	10



Sage[®] Suture Tray Packs

Suture Pack #3

PRODUCT CODE: 11011049

1 x Multipurpose Tray 220 x 120 x 20mm

5 x Non-Woven Bobs Latex Free .3gm

3 x Non-Woven Swabs 7.5 x 7.5cm 4-Ply

1 x Forceps Green

1 x Needle Holder Mayo Hegar 15cm Heavy

1 x Scissors Iris

1 x Forceps Adson Tissue 1:2

1 x Forceps Dressing Block End Blunt

1 x Disposable Scalpel #15

1 x Hydrowipe

1 x Field White Plastic 50 x 50cm

Suture Pack #6

PRODUCT CODE: 11011052

1 x Multipurpose Tray 220 x 120 x 20mm

5 x Green Gauze Swabs 10 x 10cm 8-Ply

1 x Needle Holder Mayo Hegar 15cm Heavy

1 x Scissors Sharp/Blunt

1 x Forceps Adson Tissue 1:2

1 x Forceps Dressing Block End Blunt

1 x Disposable Scalpel #11

1 x Field White Plastic 75 x 75cm

Suture Pack #7

PRODUCT CODE: 11011062

1 x Multipurpose Tray 220 x 120 x 20mm

8 x Non-Woven Swabs 7.5 x 7.5cm 4-Ply

1 x Forceps Green

1 x Needle Holder Mayo Hegar

1 x Scissors Iris

1 x Adson Tissue Forceps 1:2

1 x Forceps Dressing Block End Blunt

1 x Fenestrated Drape 4-Ply Polylined 51 x 67 cm
(Fenestration 10 x 10cm)

1 x Towel Blue Polyline 2 Ply, 42 x 40cm

1 x Field White Plastic 50 x 50cm

Suture Pack #12

PRODUCT CODE: 11011131

1 x Multipurpose Tray 220 x 120 x 20mm

5 x Non-Woven Swabs 7.5 x 7.5cm 4-Ply

3 x Gauze Swabs 10 x 10cm

1 x Forceps Green

1 x Forceps Dressing Block End Blunt

1 x Scissors Iris

1 x Forceps Adson Tissue 1:2

1 x Needle Holder Mayo Hegar 15cm Heavy

1 x Towel 2-Ply Paper 42 x 40cm

1 x Field White Plastic 50 x 50cm



Sage[®] Suture Tray Packs

Suture Pack #13

PRODUCT CODE: 11011157

1 x Multipurpose Tray 220 x 120 x 20mm
1 x Fenestrated drape 51 x 58cm (fenestration 10 x 10cm)
3 x Gauze Swabs 10 x 10cm 8-Ply
1 x Forceps Green
1 x Needle Holder Mayo Hegar 15cm (heavy)
1 x Forceps Adson Tissue 1:2
1 x Forceps Dressing Block end Blunt
1 x Scissors Iris, Closed shank
5 x Non-woven bobs Latex Free, .45gm
1 x Surgical Cotton Towel 42 x 67.5cm
1 x White Plastic Field 50 x 50cm

Suture Pack Removal

PRODUCT CODE: 11011057

1 x Standard Tray 145 x 130 x 18mm
3 x Cotton Balls .3gm
2 x Forceps Green
1 x Stitch Cutter
1 x Towel Blue Guard 42 x 42cm
1 x Field White Plastic 50 x 50cm

Suture Pack Removal #1

PRODUCT CODE: 11011133

1 x Forceps Semken Tissue
1 x Scissors Iris