



S+M[®] Under Pads Blue

S+M[®] Under Pads are made with soft, non-woven layers with a blue waterproof backing that absorbs fluids while protecting surfaces. S+M[®] Under Pads are also made with cellulose fibres which trap the fluid and with a strong, leak-proof polyethylene backing sheet. Single use ensures a hygienic and cost effective solution for incontinence care.



Suggested applications

Suitable for use in hospitals, doctors' surgeries, nursing homes and domestic households.

Suitable for use on beds, chairs, wheelchairs, vehicle seats where protective coverage is required.

Also ideal for use as a baby change mat.

Features

3 layer construction - nonwoven skin layer, 16gsm 3ply tissue layer and polyethylene waterproof backing layer
Max absorbency volume 270ml, 0.9l/min

The top non-woven layer absorbs the liquid and draws it quickly down to the second and third tissue layers, leaving the surface dryer.

The four sealed edges prevent liquids spreading beyond the surface of the pad.

Provided as bundles of 50 pieces for ease of storage and allowing single hand opening

Hygienic

Cost effective

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10126519	S+M [®] Bed Cover Under Pad 40X54cm 5ply Blue	300	50	300