

Catheter Packs

Sage Catheter Pack is carefully assembled to ensure convenience, hygiene, and ease of use during catheterisation procedures. They are designed to minimise the risk of infection and provide healthcare professionals with the essential tools needed to perform catheterisation.



Easy Open Pack

Quick and convenient one-step process



Puncture Resistant Tray

Firm structure with reinforced edges



Aids in Infection Control

Provides an on-demand sterile set-up to help minimise risk of infection



11003012



11003017

Suitable for:

- Hospitals, clinics, and other healthcare settings where catheterisation is a routine procedure

REF	Description	Sterility	Pcs/Unit	Units/Ctn
Catheter Packs				
11003001	Catheter Pack #1	Sterile	30	1
11003002	OPS™ Catheter Pack #2	Sterile	60	1
11003009	OPS™ Catheter Pack #9	Sterile	55	1
11003017	OPS™ Catheter Pack #17	Sterile	60	1
Catheter Packs with Purifil Water Syringe				
11003011	OPS™ Catheter Pack #11	Sterile	48	1
11003012	OPS™ Catheter Pack #12	Sterile	48	1
11003013	OPS™ Catheter Pack #13	Sterile	30	1
11003016	OPS™ Catheter Pack #16	Sterile	72	1
11003106	OPS™ Catheter Pack #6	Sterile	72	1
11003108	OPS™ Catheter Pack #8	Sterile	72	1

Component details

Catheter Packs

					WITH PURIFIL PREFILLED STERILE WATER SYRINGE					
REORDER	11003001	11003002	11003009	11003017	11003011	11003012	11003013	11003016	11003106	11003108
Sterile Field	75cm x 75cm White Plastic	50cm x 50cm Blue Laminated	50cm x 50cm Blue Laminated	50cm x 50cm Blue Laminated	70cm x 70cm White Plastic	70cm x 70cm White Plastic	70cm x 70cm White Plastic	70cm x 70cm White Plastic	60cm x 60cm White Plastic	60cm x 60cm White Plastic
Kidney Dish	2 x 375ml Green						1 x 700ml Clear			
OPS™ Long Tray 24.5cm x 10.5cm x 2.5cm		2 x Teal	2 x Teal	2 x Teal	2 x Teal	2 x Teal		2 x Teal	1 x Teal	1 x Teal
Gallipots 120ml	2									
Cotton Balls	20	7								5
Gauze Swabs 8ply	2 x 10cm x 10cm		5 x 7.5cm x 7.5cm		5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	6 x 10cm x 10cm	6 x 10cm x 10cm
Non-Woven Bobs				7						
Forceps	4 x Green OPS™	2 x Green OPS™	2 x Green OPS™, 1 x Yellow OPS™	2 x Green OPS™	3 x Green	3 x Green				3 x Green
Catheter Securement Device 6.5cm x 10.5cm					1		1			
Lubricating Jelly Sachet 3g					1	1	1	1	1	1
Sodium Chloride Ampoule					1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml
Prefilled Water Syringe					1 x 10ml	1 x 10ml	1 x 10ml		1 x 10ml	1 x 10ml
Nitrile Powder Free Gloves		2 x Medium		2 x Medium					2 x Large	2 x Large
Catheter Insertion Label							14cm x 10cm Yellow			
Paper Towel 40cm x 38cm	2 x 2ply	1 x 2ply	1 x 4ply	1 x 2ply					1 x 2ply	1 x 2ply
Medical Towel					30cm x 50cm	30cm x 50cm				
Fenestrated Drape 51cm x 67cm (Fen 10cm dia)		1	1	1						
Fenestrated Drape 60cm x 60cm (Fen 6cm x 10cm)									1 x Blue Cellulose/ Polyester	1 x Blue Cellulose/ Polyester
Fenestrated Drape 80cm x 100cm (Fen 9cm x 18cm)					1	1	1	1		

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

