

Catheter Packs

Sage Catheter Pack is carefully assembled to ensure convenience, hygiene, and ease of use during catheterisation procedures. They are designed to minimise the risk of infection and provide healthcare professionals with the essential tools needed to perform catheterisation with meticulous care.



Easy Open Pack

Quick and convenient one-step process



Puncture Resistant Tray

Firm structure with reinforced edges



Sterility Guaranteed

Minimising the risk of contamination



11003012



11003009

Suitable for:

- Hospitals, clinics, and other healthcare settings where catheterisation is a routine procedure

REF	Description	Sterility	Pcs/Unit	Units/Ctn
Standard Catheter Packs				
11003001	Catheter Pack #1	Sterile	30	1
11003002	Catheter Pack #2	Sterile	60	1
11003009	Catheter Pack #9	Sterile	55	1
OPS Catheter Packs with Purifil Water Syringe				
11003011	OPS™ Catheter Pack #11	Sterile	48	1
11003012	OPS™ Catheter Pack #12	Sterile	48	1
11003013	OPS™ Catheter Pack #13	Sterile	30	1
11003016	OPS™ Catheter Pack #16	Sterile	72	1
11003106	OPS™ Catheter Pack #6	Sterile	72	1
11003108	OPS™ Catheter Pack #8	Sterile	72	1

Component details

Catheter Packs

				WITH PURIFIL PREFILLED STERILE WATER SYRINGE					
REORDER	11003001	11003002	11003009	11003011	11003012	11003013	11003016	11003106	11003108
Sterile Field	75cm x 75cm White Plastic	50cm x 50cm Blue Laminated	50cm x 50cm Blue Laminated	70cm x 70cm White Plastic	70cm x 70cm White Plastic	70cm x 70cm White Plastic	70cm x 70cm White Plastic	60cm x 60cm White Plastic	60cm x 60cm White Plastic
Kidney Dish	2 x 375ml Green					1 x 700ml Clear			
OPSTM Long Tray 24.5cm x 10.5cm x 2.5cm		2 x Teal	2 x Teal	2 x Teal	2 x Teal		2 x Teal	1 x Teal	1 x Teal
Gallipots 120ml	2								
Cotton Balls	20	7							5
Gauze Swabs 8ply	2 x 10cm x 10cm		5 x 7.5cm x 7.5cm	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	5 x 7.5cm x 7.5cm XRD	6 x 10cm x 10cm	6 x 10cm x 10cm
Forceps	4 x Green OPSTM	2 x Green OPSTM	2 x Green OPSTM, 1 x Yellow OPSTM	3 x Green	3 x Green				3 x Green
Catheter Securement Device 6.5cm x 10.5cm				1		1			
Lubricating Jelly Sachet 3g				1	1	1	1	1	1
Sodium Chloride A mpoule				1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml	1 x 30ml
Prefilled Water Syringe				1 x 10ml	1 x 10ml	1 x 10ml		1 x 10ml	1 x 10ml
Nitrile Powder Free Gloves		2 x Medium						2 x Large	2 x Large
Catheter Insertion Label						14cm x 10cm Yellow			
Paper Towel 40cm x 38cm	2 x 2ply	1 x 2ply	1 x 4ply					1 x 2ply	1 x 2ply
Medical Towel				30cm x 50cm	30cm x 50cm				
Fenestrated Drape 51cm x 67cm (Fen 10cm x 10cm)		1	1						
Fenestrated Drape 60cm x 60cm (Fen 6cm x 10cm)								1 x Blue Cellulose/ Polyester	1 x Blue Cellulose/ Polyester
Fenestrated Drape 80cm x 100cm (Fen 9cm x 18cm)				1	1	1	1		

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

