

## Catheter Securement Devices

Sage Catheter Securement Devices utilise a combination of adhesive and securement straps to effectively secure the catheter in place. The adhesive backing provides a secure and durable bond to the skin, while the straps offer additional support and stability.



### Re-adjustable Tabs

Ensures proper positioning and stability of catheters



### Water Resistant

Allowing patients to maintain their regular hygiene routines without compromising the device's effectiveness




### Gentle Adhesion


Easy application and removal, minimising discomfort for patients

**Suitable for:**


- Urinary catheter drainage
- Universal securement for various tubing and drains




11020001



11020002



11020003



11020004



REF	Size	Sterility	Pcs/Unit	Units/Ctn
11020001	10.5cm x 6.5cm	Sterile	50	4
11020002	16.8cm x 5.7cm	Sterile	50	4
11020003	9cm x 3.5cm	Sterile	50	4
11020004	12.5cm x 8.5cm	Sterile	50	4

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



**Sentry Medical Pty Ltd**  
 1300 995 999 | sales@sentrymedical.com.au  
 www.sentrymedical.com.au

