

IV Starter Pack #9

Convenient and easy to use, the Sage IV Starter Pack #9 provides an all-in-one solution with essential items to initiate intravenous therapy swiftly and effectively. This pack includes Asguard IV Dressing Reinforced 7cm x 8.5cm (ISB303). The IV Starter Pack #9 belos maintain a sterile environment, ensuring that healthcare professionals can perform the

The IV Starter Pack #9 helps maintain a sterile environment, ensuring that healthcare professionals can perform the procedure safely and hygienically, reducing the risk of infections.



All-in-one Pack

Comprehensive components help streamline medical procedures



Aids in Infection Control

Provides an on-demand sterile set-up to help minimise risk of infection



Clinical and Time Efficient

Access to vital components for IV therapy and easy setup





Comes in a convenient desktop dispenser of 45 pieces.

Suitable

for:

- Outpatient clinics
- Surgical procedures
- Emergency rooms and trauma care

REF	Description	Components	Sterility	Pcs/Unit	Units/Ctn
11005042	IV Starter Pack #9	1 x Blue Polylined Towel - 42cm x 40cm 1 x IV Starter Pack Tray 1 x Swab Stick, 2% CHG, 70% Isopropyl Alcohol 3 x Non-Woven Swabs - 7.5cm x 7.5cm, 4ply 1 x Asguard IV Dressing Reinforced 7cm x 8.5cm (ISB303) 2 x Patient Label - Green Relocated IV Site 1 x White Plastic Field 47.5cm x 47.5cm	Sterile	180	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.









Asguard[®]



IV Dressing Reinforced

The Asguard Reinforced IV Dressings are constructed from a High MVTR film, which allows for excess moisture to evaporate, combined with a laminated non-woven border for securement and enhanced fixation. They maintain firm adhesion and allow for insertion site monitoring.



High MVTR

Provides environment with optimal moisture management



Secure Adhesion

Adhesive dressing enables firm fixation, minimising device movement or premature displacement



Waterproof

Protective layer prevents liquid from passing through covering







Suitable

- Protection and securement of PIVC lines

fo	r:	

REF	Size	Туре	Sterility	Pcs/Unit	Units/Ctn
ISB303	7cm x 8.5cm	Bordered	Sterile	100	18

Component details



Asguard°



IV Film Dressing Range



REORDER	SIZE	SUITABLE FOR	STERILITY	PCS/UNIT	UNITS/CTN
ISB301	5cm x 5.7cm	Paediatric, PIVC	Sterile	100	10
ISB302	6.5cm x 7cm	PIVC	Sterile	100	20
ISB303	7cm x 8.5cm	PIVC	Sterile	100	18
ISB304	7cm x 9cm	PIVC	Sterile	100	12
ISB312	8.5cm x 10.5cm	PIVC, CVC	Sterile	100	16
ISB402	6cm x 7cm	Paediatric, PIVC	Sterile	100	24
ISB403	10cm x 12cm	PICC, PORTS, PIVC, CVC	Sterile	100	12
ISB404	11cm x 15cm	PICC, PORTS, PIVC, CVC	Sterile	100	12
ISB202	6cm x 7cm	PICC, PORTS, PIVC, CVC	Sterile	50	10
ISB203	10cm x 12cm	PICC, PORTS, PIVC, CVC	Sterile	25	10
IVD002	6cm x 7cm	PICC	Sterile	50	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.







