

## Dust Cover Bags

Sentry Dust Cover Bags are ideal for storage of sterilised instruments and trays, protecting them from dust during storage or transportation. For extra protection, the bags can be heat sealed using an impulse bar sealer. For sealing, please follow the heat sealer manufacturer’s instructions and facility protocol.



**Heat Sealable Bags**  
Easy to seal using an impulse bar sealer



**Transparent**  
Clear for ease of viewing contents



**Protection**  
Protects sterilised instruments and trays from dust



- Suitable for:**
- Storage of sterilised instruments and trays
  - Extending sterile tenability

REF	Size	Colour	Pcs/Unit	Units/Ctn
DC001	152mm x 254mm	Clear	100	20
DC002	203mm x 305mm	Clear	100	15
DC003	254mm x 381mm	Clear	100	15
DC004	305mm 508mm	Clear	100	15
DC005	406mm x 508mm	Clear	100	5
DC006	450mm x 762mm	Clear	100	4
DC007	610mm x 762mm	Clear	100	4

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

