




SureSafe CliniCover

CliniCover is used in the operating room to protect the operating table where a fully absorbent underlay is required. CliniCover can replace regular linen and it has an impervious backing preventing blood and body fluids from coming into contact with the operating table and being absorbed by mattress pad.



Features

-
- Lint free
-
- Latex free
-
- Strong absorbent reinforcement
-
- Impervious backing
-
- Disposable
-

Benefits

-
- Reduced risk of contamination
-
- Improved infection control
-
- Easy of use
-
- Reduced cleaning time
-
- Greater patient comfort
-
- Extends the life of mattress and operating table pads
-

FABRIC KEY	 Absorbent Impervious Fabric	 Reinforcement (strong absorbent Fabric)	 Incise Film	 Side Plastic	 Clear Plastic Film	 Instrument Pad	 Adhesive Area
-------------------	---------------------------------	---	-----------------	------------------	------------------------	--------------------	-------------------

PRODUCT CODE	PRODUCT DESCRIPTION	SIZE	SKU/TU (MOQ)	SKU/CTN
S60040	Medium, LW	60 x 90cm	120	120
S60010	Large	100 x 225cm	40	40
S60030	Large, LW	120 x 228cm	40	40