



SureSafe Orthopaedic Hand Drape

Hand Drape 196 x 371cm

SureSafe 35 draping material


Strong absorbent reinforcement 42 x 85cm

5cm Elastic fenestration

Gusseted hand table extension

Velcro tube and cord holders

Sterile individually packed



SureSafe Draping Material

3-layer polypropylene for anti-static, fluid and abrasion resistance including 2 outer layers spunbond polypropylene and 1 inner layer of meltblown polypropylene

Low-linting with a laminated non-skid surface








Strong and safe fluid control features, soaks up spill and splatters

No fluid strike through, fully impervious drapes

Absorbent throughout the whole drape

Easy handling and minimal size footprint folding techniques

Compliant and tested to EN 13795

FABRIC KEY							
	Absorbent Impervious Fabric	Reinforcement (strong absorbent Fabric)	Incise Film	Side Plastic	Clear Plastic Film	Instrument Pad	Adhesive Area


PRODUCT CODE	PRODUCT DESCRIPTION	OUTER DRAPE DIMENSION (CM)	REINFORCED DRAPE DIMENSION (CM)	SKU/TU (MOQ)	SKU/CTN
S72603	Hand Drape	196 x 371	42 x 85	8	8



SureSafe Orthopaedic Shoulder Drape

Beach Chair Shoulder Drape

-
- Total Shoulder Replacement draping solution
-
- Strong absorbent reinforcement 53 x 66cm
-
- Tear drop fenestration
-
- Attached fluid contract pouch with exit port
-
- Velcro tube and cord holders
-
- Sterile individually packed
-



S72626

SureSafe Draping Material

-
- 3-layer polypropylene for anti-static, fluid and abrasion resistance including 2 outer layers spunbond polypropylene and 1 inner layer of meltblown polypropylene
-
- Low-linting with a laminated non-skid surface
-
- Strong and safe fluid control features, soaks up spill and splatters
-
- No fluid strike through, fully impervious drapes
-
- Absorbent throughout the whole drape
-
- Easy handling and minimal size footprint folding techniques
-
- Compliant and tested to EN 13795
-

FABRIC KEY	 Absorbent Impervious Fabric	 Reinforcement (strong absorbent Fabric)	 Incise Film	 Side Plastic	 Clear Plastic Film	 Instrument Pad	 Adhesive Area
------------	---------------------------------	---	-----------------	------------------	------------------------	--------------------	-------------------

PRODUCT CODE	PRODUCT DESCRIPTION	OUTER DRAPE DIMENSION (CM)	REINFORCED DRAPE DIMENSION (CM)	SKU/TU (MOQ)	SKU/ CTN
S72626	Beach Chair Shoulder Drape	259 x 406	53 x 66	10	10





SureSafe Orthopaedic Head Drape Pack

Head Drapes Pack

- 1 ENT split drape 193 x 315cm
- 2 impervious drapes 70 x 45cm heat joined
- 1 standard back table cover 112 x 228cm
- Strong absorbent reinforcement 50 x 45cm
- Instrument pad
- Velcro tab holders
- Sterile individually packed

SureSafe Draping Material

- 3-layer polypropylene for anti-static, fluid and abrasion resistance including 2 outer layers spunbond polypropylene and 1 inner layer of meltblown polypropylene
- Low-linting with a laminated non-skid surface
- No fluid strike through, fully impervious drapes
- Absorbent throughout the whole drape
- Easy handling and minimal size footprint folding techniques
- Compliant and tested to EN 13795



PRODUCT CODE	PRODUCT DESCRIPTION	ENTT SPLIT SIZE (CM)	REINFORCED AREA DIMENSION (CM)	SKU/TU (MOQ)	SKU/ CRT
SC82805	Head Drapes Pack	7.6 x 63.5	50 x 45	12	12