

# SureSafe Scrub Shirts and Pants

## Scrub Shirts

3 layer SMS material certified AAMI Level 2

Suitable for light fluid applications

Soft, comfortable and breathable

Round neck, reinforced pocket

Non-sterile, available in three sizes



## Scrub Pants


3 layer SMS material certified AAMI Level 2

Suitable for light fluid applications

Soft, comfortable and breathable

Reinforced seams, durable drawstrings

Non-sterile, available in three sizes



PRODUCT CODE	PRODUCT DESCRIPTION	SIZE	COLOUR	SKU/TU (MOQ)	SKU/CTN
PA51001M	SureSafe Scrub Shirt	Medium	Blue	50	50
PA51001L	SureSafe Scrub Shirt	Large	Blue	50	50
PA51001XL	SureSafe Scrub Shirt	X-Large	Blue	50	50
PA52001M	SureSafe Scrub Pant	Medium	Blue	50	50
PA52001L	SureSafe Scrub Pant	Large	Blue	50	50
PA52001XL	SureSafe Scrub Pant	X-Large	Blue	50	50