



SureSafe Standard Surgical Gown

SureSafe® standard surgical gowns are made from 4-layer polypropylene materials with anti-static and fluid resistance. Two layers of spunbond polypropylene ensure an extraordinary flexibility and drapeability, provide tough abrasion resistance and cloth like feel, and two layers of meltblown polypropylene forming an exceptional microbial barrier against bacteria.

Compliance AAMI Level 2 - EN 13795 – Full Coverage

4-layer anti-static breathable fabric

Fluid repellent

Excellent microbial barrier

Light weight material, less recycled waste

Velcro neck binding

Cost-effective

Generous cut

Large transfer wrap – 90 x 90cm

2 paper towels standard or 1 cotton huck towel

Available non-sterile

Sterile individually packed, available non-sterile*




*Contact Aaxis sales consultant for lead time and MOQ

PRODUCT CODE	MATERIAL	REINFORCED	SIZE	PAPER / COTTON TOWEL	AAMI LEVEL	SKU/TU (MOQ)	SKU/ CARTON
S71652P	35gsm 4 Layer SMMS	No	Medium	2 paper	Level 2	36	36
S71653P	35gsm 4 Layer SMMS	No	Large	2 paper	Level 2	36	36
S71654P	35gsm 4 Layer SMMS	No	X-Large	2 paper	Level 2	36	36
S71655P	35gsm 4 Layer SMMS	No	XX-Large	2 paper	Level 2	36	36

SureSafe Standard Surgical Gown

Product Dimensions



PRODUCT CODE	SIZE	DIMENSIONS								
		A	B	C	D	E	F	G	H	I
S71652	Medium	120	72	112	55	26	22	57	8	150
S71653	Large	128	72	118	58	27	24	65	8	156
S71654	Extra Large	143	72	134	60	28	24	64	8	160
S71655	Extra-Extra Large	148	85	140	65	32	26	77	8	175
	Tolerance	±2.0	±1.5	±2.0	±1.5	±1.0	±1.0	±1.5	±0.75	±2.0