



SureSafe Standard Surgical Gown

SureSafe® standard surgical gowns are made from 5-layer polypropylene materials with anti-static and fluid resistance. Three layers of spunbond polypropylene ensure an extraordinary flexibility and drapeability, provide tough abrasion resistance and cloth like feel, and two layers of meltblown polypropylene forming an exceptional microbial barrier against bacteria.

5-layer anti-static breathable fabric - SMMMS

Compliance AAMI Level 3 - EN 13795 – Full Coverage

Fluid repellent

Excellent microbial barrier

Excellent physical barrier to contaminants/fluids

Stronger fabric, matte finish, resistance to tears

Softer fabric with comfort inside and outside

Generous cut

Velcro neck binding

Book fold

Sterile individually packed, available non-sterile*




*Contact Aaxis sales consultant for lead time and MOQ

PRODUCT CODE	MATERIAL	REINFORCED	SIZE	PAPER / COTTON TOWEL	AAMI LEVEL	SKU/TU (MOQ)	SKU/ CARTON
S71152	50gsm 5 Layer SMMMS	No	Medium	1 cotton	Level 3	24	24
S71153	50gsm 5 Layer SMMMS	No	Large	1 cotton	Level 3	24	24
S71154	50gsm 5 Layer SMMMS	No	XL-Large	1 cotton	Level 3	24	24
S71155	50gsm 5 Layer SMMMS	No	XXL-Large	1 cotton	Level 3	24	24

SureSafe Standard Surgical Gown

Product Dimensions



PRODUCT CODE	SIZE	DIMENSIONS								
		A	B	C	D	E	F	G	H	I
S71152	Medium	120	72	112	55	26	22	57	8	150
S71153	Large	128	72	118	58	27	24	65	8	156
S71154	Extra Large	143	72	134	60	28	24	64	8	160
S71155	Extra-Extra Large	148	85	140	65	32	26	77	8	175
	Tolerance	±2.0	±1.5	±2.0	±1.5	±1.0	±1.0	±1.5	±0.7	±2.0