

Fenestrated Drape

Suresafe Fenestrated Drapes create an absorbent and impervious barrier to fluid strikethrough. They are designed to provide a protective covering to isolate a surgical incision or surgical site. These drapes come in a variety of fenestration sizes and shapes.



Strong and Durable

Made from low linting and abrasion resistant material





Absorbent

Non-woven absorbent outer layer to maintain a dry working area



Impervious

Helps protect against fluid and contaminant strikethrough to aid in infection control



Suitable for:

- Operating theatres
- Emergency departments
- Procedure packs

| REF | Size | Hole Size | Hole Shape | Sterility | Pcs/Unit | Units/Ctn |
|-----------------|---------------|--------------|------------|-----------|----------|-----------|
| S72221 | 45cm x 60cm | 6cm x 7.5cm | Oval | Sterile | 120 | 1 |
| S72222 | 45cm x 75cm | 6cm x 7.5cm | Oval | Sterile | 120 | 1 |
| S72223 | 60cm x 60cm | 6cm x 7.5cm | Oval | Sterile | 100 | 1 |
| S72226 | 90cm x 90cm | 6cm x 7.5cm | Oval | Sterile | 100 | 1 |
| S72230 | 112cm x 150cm | 10cm x 15cm | Oval | Sterile | 20 | 1 |
| S72232 | 170cm x 270cm | 25cm x 25cm | Oval | Sterile | 30 | 1 |
| S72291 | 60cm x 60cm | 2.5cm x 6cm | Rectangle | Sterile | 200 | 1 |
| S72292 | 60cm x 60cm | 3.5cm x 10cm | Rectangle | Sterile | 160 | 1 |
| S72293 | 60cm x 90cm | 3.7cm x 11cm | Rectangle | Sterile | 160 | 1 |
| S72294 | 90cm x 122cm | 6cm x 15cm | Rectangle | Sterile | 100 | 1 |
| 11010012 | 50cm x 68cm | 10cm x 10cm | Round | Sterile | 250 | 1 |

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

