

## Sterilisation Wraps

Suresafe Sterilisation Wraps are made from non-woven, multi-layer SMS material and are intended for wrapping surgical trays and instruments before they undergo the process of sterilisation. The Suresafe Sterilisation Wraps provide a physical and microbial barrier that helps maintain sterility.



### Range of Sizes

Offers versatility to the dynamic demands of the CSSD environment




### Low Linting

Aids in infection control and minimises the risk of contamination



### Suitable for Most Sterilisation Methods

Can be used in steam, ethylene oxide, or plasma sterilisation processes



### Suitable for:

- Light to medium tray
- Linen packs
- Medical devices & instruments
- Basin sets

REF	Size	Weight	Sterility	Pcs/Unit	Units/Ctn
<b>J6102</b>	45cm x 45cm	45gsm	Non-Sterile	500	1
<b>J6103</b>	60cm x 60cm	45gsm	Non-Sterile	500	1
<b>J6104</b>	76cm x 76cm	45gsm	Non-Sterile	300	1
<b>J6105</b>	91cm x 91cm	45gsm	Non-Sterile	150	1
<b>J6106</b>	101cm x 101cm	45gsm	Non-Sterile	150	1
<b>J6107</b>	94cm x 140cm	45gsm	Non-Sterile	100	1
<b>J6108</b>	121cm x 121cm	45gsm	Non-Sterile	100	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



**Sentry Medical Pty Ltd**

1300 995 999 | sales@sentrymedical.com.au  
www.sentrymedical.com.au

