

## Sterilisation Wraps

Suresafe Sterilisation Wraps are crafted from 100% polypropylene with a robust 5 layer SMMMS construction ensuring optimal protection for medical instruments and equipment.  
Available in a range of weights and sizes offering versatility to the dynamic demands of the CSSD environment.



### Australian Made

Locally manufactured In Australia




### Low Linting

Aids in infection control and minimises the risk of contamination



### Suitable for Moist Sterilisation Methods

Can be used in steam, ethylene oxide, or plasma sterilisation processes



### Suitable for:

#### Standard

- Light to medium tray
- Linen packs
- Medical devices & instruments
- Basin sets

#### Heavy Duty

- Large trays
- Instrument sets
- Medical devices
- Multiple basin sets

REF	Description	Size	Weight	Pcs/Unit	Units/Ctn
<b>STANDARD</b>					
J6102	Suresafe Wrap - Standard	45cm x 45cm	47gsm	500	1
J6103	Suresafe Wrap - Standard	60cm x 60cm	47gsm	500	1
J6104	Suresafe Wrap - Standard	76cm x 76cm	47gsm	300	1
J6105	Suresafe Wrap - Standard	91cm x 91cm	47gsm	150	1
J6106	Suresafe Wrap - Standard	101cm x 101cm	47gsm	150	1
J6107	Suresafe Wrap - Standard	94cm x 140cm	47gsm	100	1
J6108	Suresafe Wrap - Standard	121cm x 121cm	47gsm	100	1
<b>HEAVY DUTY</b>					
J6201	Suresafe Wrap - Heavy Duty	76cm x 76cm	63gsm	200	1
J6202	Suresafe Wrap - Heavy Duty	91cm x 91cm	63gsm	150	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

