

ThreadX®

X-ray Detectable Abdominal Sponge

ThreadX®



Low
Linting



6 Ply



High Fluid
Absorbency

ThreadX® Abdominal Sponges are without tapes and constructed from 100% cotton layered gauze, making them highly absorbent and non abrasive. A black x-ray detectable thread is internally heat bonded to the gauze.

HIGHLY ABSORBENT

Approximate absorption capacity over 8 times its own weight in fluid

HEMMED SEAMS

Preventing fraying of edges

X-RAY DETECTABLE

X-ray detectable thread for patient safety



Suitable for:

All open surgical procedures

Operating room

Reorder	Ply	Size	Pcs	Sterility	Pcs/Unit	Units/Ctn
ABX002	6ply	22.5cm x 22.5cm	5pcs	Sterile	60	1
ABX004	6ply	30cm x 30cm	5pcs	Sterile	50	1
ABX008	6ply	50cm x 30cm	5pcs	Sterile	30	1

