

ThreadX[®]

Hybrid X-ray Detectable Abdominal Sponge

ThreadX[®]



Low
Linting



5ply



Very High Fluid
Absorbency

ThreadX Hybrid XRD Abdominal Sponges are without tapes and constructed from non-woven inner layers with an all-cotton gauze outer cover. A black x-ray detectable thread is heat bonded to the gauze.

X-RAY DETECTABLE

Black x-ray detectable thread for patient safety

HIGHLY ABSORBENT

Approximate absorption capacity over 10 times its own weight in fluid

HEMMED SEAMS

Preventing fraying to edges



Suitable for:

Open surgical procedures where improved fluid control is advantageous

Reorder	Ply	Size	Pcs	Sterility	Pcs/Unit	Units/Ctn
ABX112	5ply	22.5cm x 22.5cm	5pcs	Sterile	60	1
ABX114	5ply	30cm x 30cm	5pcs	Sterile	50	1
ABX103	5ply	30cm x 50cm	5pcs	Non Sterile	32	1
ABX118	5ply	30cm x 50cm	5pcs	Sterile	30	1



Sterile Packaging