

# Urimaax<sup>®</sup> A1 Drainage Bags



10013200



10013101

## Features

Latex free

Non sterile

Push-pull bottom drain valve

Anti-reflux valve incorporated into bag

Universal stepped catheter connector

120cm kink resistant tubing with regular bore

Clear, graduated front with white frosted backing to allow easy visual readings

Hanger holes fit universal plastic or metal hangers

100ml graduated measurements for easy identification of volume

2000ml capacity

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10013101	Urimaax <sup>®</sup> Drainage Bag B/Drain Push-Pull A1 2000mls	200	20	200
10013200	Urimaax <sup>®</sup> Drainage Bag Accessories Bag Hanger Plastic	100	10	100