

Urimaax® A3 & A3+ Standard Drainage Bags



10013120



10013200



10013103

Features

Latex free

Non sterile

Anti-reflux valve incorporated into bag

Universal stepped catheter connector

Clear, graduated front with white frosted backing to allow easy visual readings

Hanger holes fit universal plastic or metal hangers

100ml graduated measurements for easy identification of volume

2000ml capacity

Blue dual plastic hanger available

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10013120	Urimaax® Drainage Bag B/Drain T-Bar A3 Standard 2000mls	200	10	200
10013103	Urimaax® Drainage Bag B/Drain T-Bar A3 2000mls	200	10	200
10013200	Urimaax® Drainage Bag Accessories Bag Hanger Plastic	100	10	100