

# Urimaax® Drainage Bags Leg

## Features

- 
- Kink resistant tubing
- 
- Tubing from 60mm to 500mm
- 
- Latex-free velcro straps and bag
- 
- Fast and safe flip-tap bottom drain
- 
- Anti-reflux valve incorporated into bag
- 
- Soft cloth backing for added comfort
- 
- Universal stepped catheter connector
- 
- 100ml graduated measurements for easy identification of volume
- 
- Available in 500ml and 750 ml capacities
- 



## STERILE

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10013302	Urimaax®*** Drainage Bag Leg 500ml 6cm Sterile	80	20	80
10013303	Urimaax®*** Drainage Bag Leg 500ml 10cm Sterile	80	20	80
10013304	Urimaax®*** Drainage Bag Leg 500ml 30cm Sterile	80	20	80
10013305	Urimaax®*** Drainage Bag Leg 500ml 50cm Sterile	80	20	80
10013402	Urimaax®*** Drainage Bag Leg 750ml 6cm Sterile	80	20	80
10013403	Urimaax®*** Drainage Bag Leg 750ml 10cm Sterile	80	20	80
10013404	Urimaax®*** Drainage Bag Leg 750ml 30cm Sterile	80	20	80
10013405	Urimaax®*** Drainage Bag Leg 750ml 50cm Sterile	80	20	80

## NON - STERILE

PRODUCT CODE	PRODUCT DESCRIPTION	SKU/TU (MOQ)	SKU/ Inner	SKU/ CTN
10013503	Urimaax®*** Drainage Bag Leg 500ml 10cm	80	20	80
10013504	Urimaax®*** Drainage Bag Leg 500ml 30cm	80	20	80
10013603	Urimaax®*** Drainage Bag Leg 750ml 10cm	80	20	80
10013604	Urimaax®*** Drainage Bag Leg 750ml 30cm	80	20	80
10013110	Urimaax®*** Drainage Bag Night 2000mls	250	10	250