



Easy To Use



Disposable



Infection Control

Sentry Medical Emesis Bags have a firm polypropylene twist and lock rim to provide easy control, and help prevent spillage after use. They isolate splash back and have a graduated indicative volume system to 1500ml, allowing measurement of contents in the bag. On completion of use, the bag is easily disposed of with no mess. Simply twist the bag and lock into the keyhole slot on the rim.

RED RIMS

Rims are colour coded red for infection control

INFECTION CONTROL

Easy non-leak twist and lock seal for a no mess disposal

CONVENIENT DISPENSERS

Handy wall dispenser available for hospital setting, or 3 pack dispenser for retail environments



Suitable for:

Infection controlled area

First aid

Hospitals and medical

Reorder	Description	Sterility	Pcs/Unit	Units/Ctn
VB002	VBag® Emesis Bag (Red Ring)	Non-Sterile	50	10
VB011	Sickeze® Emesis Bag Retail Shelf Pack (3pcs/pk - 6pks)	Non-Sterile	10	1