

Pillow Sleeves

Verasoft Pillow Sleeves are disposable, non-woven covers made from spunbond polypropylene fabric. The covers are intended to be used as protective covers for pillows to aid in infection control.



Non-Abrasive Surface
Soft and comfortable for patient skin



Disposable
Aids in maintaining infection control and reducing excess laundering



Ideal for Infection Control
Protects pillows from contaminants including skin cells, hair, and bodily fluids



Suitable for:

- Hospital and medical
- Infection control

REF	Size	Colour	Sterility	Pcs/Unit	Units/Ctn
10125511	70cm x 45cm	White	Non-Sterile	200	1
10125512	70cm x 45cm	Navy Blue	Non-Sterile	200	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

