



Conforms
to Body



Comfortable



Cost
Effective

The Medicrepe (Elastic Crepe Bandage) Heavy range is constructed from a blend of cotton and polyurethane fibres (spandex), giving excellent stretch and memory. Ideal for providing support and stability whilst still allowing movement.

STRONG SUPPORT

Helps immobilise body part to help prevent further injury or pain

FAST EDGES

Prevents bandage from fraying and losing shape

EASY TO USE

Simply hold the bandage in one hand and wrap circumferentially from the outside inward



Suitable for:

Strong support bandaging

Joint injuries

First Aid

Ligament injuries and fractures

Reorder	Size	Feature	Sterility	Pcs/Unit	Units/Ctn
WB002H	5cm x 1.5m	Unstretched	Non-Sterile	12	60
WB007H	5cm x 2.3m	Unstretched	Non-Sterile	12	60
WB003H	7.5cm x 1.5m	Unstretched	Non-Sterile	12	40
WB008H	7.5cm x 2.3m	Unstretched	Non-Sterile	12	40
WB004H	10cm x 1.5m	Unstretched	Non-Sterile	12	30
WB009H	10cm x 2.3m	Unstretched	Non-Sterile	12	30
WB005H	15cm x 1.5m	Unstretched	Non-Sterile	12	20

