



Highly Absorbent



100% Cotton



X-Ray Detectable

ThreadX Throat Packs are constructed from 100% cotton. A green x-ray detectable thread is heat bonded to the gauze with a black tail.

HIGHLY ABSORBENT

Suitable for controlling saliva and bleeding

X-RAY DETECTABLE WITH TAIL

For ease of visibility and patient safety



Suitable for:

Open surgical procedures involving the throat region

Prevent and/or control bleeding

Reorder	Ply	Size	Sterility	Pcs/Unit	Units/Ctn
XRT001	4ply	5cm x 1m	Sterile	100	1
XRT002	4ply	5cm x 2m	Sterile	100	1

