

Sentry Medical Zinc Oxide Tape



Durable



Contours
To Body



Water
Resistant

Sentry Medical Zinc Oxide Tapes are constructed from plain cotton weave fabric with zinc oxide adhesive. This tape assists in the prevention of soft tissue damage as well as supporting minor sprains or strains.

RIGID AND STRONG

Aids in keeping bandages in place

ZINC OXIDE ADHESIVE

Helps adhere to moist or sweaty skin

Suitable for:

Securing dressings

Sport injuries

Protecting Skin (e.g. blisters)

Sprains or strains



Reorder	Size	Pcs/Unit	Units/Ctn
ZAT001	2.5cm x 5m	12	30
ZAT002	1.25cm x 5m	24	30