

Adson Forceps - Sterile

Sage Adson Forceps are crafted from high-quality stainless steel, ensuring durability and longevity. The fine tips provide a secure grip, allowing for precise manipulation of delicate tissues and structures during surgical procedures.



Secure Grip

Allowing for a secure and controlled grip



Fine Tips

For precise grasping and manipulation of tissues



Blue Handles

Easy to identify as single use



Toothed



No Teeth

Suitable for:

- Dissecting, clamping, and holding tissue in place during suturing
- Grasping and manipulating tissues during medical procedures

REF	Size	Feature	Sterility	Pcs/Unit	Units/Ctn
11011018	12cm	No Teeth	Sterile	60	1
11011019	12cm	1:2 Teeth	Sterile	60	1
11011165	12cm	No Teeth, Serrated	Sterile	50	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

