


Mosquito Halsted Forceps - Sterile

Sage Mosquito Halsted Forceps are designed with a slender, elongated shape and fine jaws that taper to a point. The forceps are made of high-quality stainless steel, ensuring durability and corrosion resistance.



Finger-ring Handles

Comfortable grip allowing for precise control



Fine Tips

For precise and controlled handling of small objects or tissues



Blue Handle

Easy to identify as single use

Suitable for: - Controlling bleeding in surgical procedures



REF	Size	Feature	Sterility	Pcs/Unit	Units/Ctn
11011020	12.5cm	Straight	Sterile	25	1
11011022	12.5cm	Curved	Sterile	20	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

