

Flaminal® Enzyme Alginogel

Flaminal® is an enzyme alginogel – a wound healing agent that contains a hydrated alginate with an embedded antibacterial enzyme system. It can be used for acute and chronic wound types, or wounds at risk of infection: Flaminal® Hydro for slightly to moderately exuding wounds and Flaminal® Forte for wounds with moderate to high amounts of wound exudate.



Debridement Gel

Dissolves necrotic tissue and fibrinous slough and absorbs them from the wound bed



Absorbent Alginate

Absorbs excess exudate, enables autolytic debridement (cleaning of wound and absorption of bacteria)



Antimicrobial Enzyme System

Offers antimicrobial protection, and reduces bacteria released from a biofilm



Suitable for:

- Superficial and deep partial thickness burns
- Diabetic foot ulcers, leg ulcers or pressure ulcers
- Infected wounds
- Neonates/paediatrics, oncologic, post-surgical or traumatic wounds
- Skin tears

REF	Description	Sterility	Pcs/Unit	Units/Ctn
1015	Flaminal Forte 50g	Non-Sterile	1	120
1016	Flaminal Forte 500g	Non-Sterile	1	18
1021	Flaminal Hydro 50g	Non-Sterile	1	120
1010	Flaminal Hydro 500g	Non-Sterile	1	18

Flaminal® should only be used under the supervision of a healthcare professional. If you are concerned about your wound, consult a healthcare professional.



Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

