

Asguard[®]

HI-TAC DRESSING RANGE



NEW

Choosing the Right IV Dressing is Important

For over 30 years, Sentry Medical has been helping health care professionals provide exceptional care to their patients through a modern and quality driven portfolio.

In today's healthcare landscape, the selection of intravenous (IV) dressings is more than a routine task—it's a critical decision that directly impacts patient safety, clinical efficiency, and healthcare costs.



DRESSING RANGE OVERVIEW



ISB301HT



ISB302HT



ISB303HT



ISB304HT



ISB312HT



ISB402HT



ISB403HT/4HT



ISB202HT/3HT



**Infection
Prevention**



**Securement
& Stability**



**Skin Integrity
& Patient Comfort**



**Aligned with
Clinical Guidelines**

HI-TAC IV Dressings:

Secure. Transparent. Trusted.

The **High-Adhesion Transparent IV Dressings - Asguard IV Hi-Tac** is designed to meet the highest standards of patient care and clinical efficiency.

In fast-paced clinical environments, securement is critical. High-adhesion IV dressings are designed to provide reliable fixation, reduce complications, and support best practices in vascular access management.



Strong Securement Reduces Dislodgement

High-adhesion dressings maintain stability in high-movement areas, lowering reinsertion and infection risks



Extended Wear Time Improves Efficiency

Designed for up to 7 days of wear, reducing dressing changes and saving nursing time



Skin-Friendly Adhesion

Modern adhesives balance strong hold with atraumatic removal to protect patients with fragile or compromised skin



**SCAN HERE
FOR PRODUCT
INFORMATION**

HI-TAC Transparent Film Dressings

Secure Protection with Clear Visibility.

The **High-Adhesion Transparent Film Dressings - Asguard Hi-Tac** is designed to meet the highest standards of patient care and clinical efficiency.

The Asguard Film Dressing HT is a transparent polyurethane adhesive film with moisture vapour permeability which has Hi-Tac adhesion to the skin and forms a flexible and comfortable layer.



Hi-Tac Secure Adhesion

Adhesive dressing enables firm fixation, minimising dressing movement or premature displacement



Waterproof

Outer layer protects site from fluid penetration, allowing patient to shower with the dressing in-situ



Transparent Film


Allows visualisation of wound, incision or IV Site



**SCAN HERE
FOR PRODUCT
INFORMATION**



REF	DESCRIPTION	SIZE	PCS/UNIT	UNITS/CTN
ISB001HT	Asguard Transparent Film Dressing Hi-Tac	4cm x 5cm	50	10
ISB002HT	Asguard Transparent Film Dressing Hi-Tac	6cm x 7cm	50	10
ISB003HT	Asguard Transparent Film Dressing Hi-Tac	10cm x 12cm	25	10
ISB004HT	Asguard Transparent Film Dressing Hi-Tac	15cm x 20cm	25	10



REF	SIZE	SUITABLE FOR	PCS/UNIT	UNITS/CTN
HI-TAC IV REINFORCED DRESSING				
ISB301HT	5cm x 5.7cm	Paediatric, PIVC	100	10
ISB302HT	6.5cm x 7cm	PIVC	100	20
ISB303HT	7cm x 8.5cm	PIVC	100	18
ISB304HT	7cm x 9cm	PIVC	100	12
ISB312HT	8.5cm x 10.5cm	PIVC, CVC	100	16
ISB402HT	6cm x 7cm	Paediatric, PIVC	100	24
ISB403HT	10cm x 12cm	PICC, PORTS, PIVC, CVC	100	12
ISB404HT	11cm x 15cm	PICC, PORTS, PIVC, CVC	100	12
HI-TAC IV DRESSING				
ISB202HT	6cm x 7cm	PICC, PORTS, PIVC, CVC	50	10
ISB203HT	10cm x 12cm	PICC, PORTS, PIVC, CVC	25	10