

IV Dressing Reinforced HT

Asguard IV Dressing Reinforced HT are constructed from transparent polyurethane adhesive film which has Hi-Tac adhesion that enables firm adhesion for secure fixation. Selected dressings incorporate a laminated non-woven reinforcement border to enhance securement while the clear film allows monitoring of the insertion site without removing the dressing.



Hi-Tac Secure Adhesion

Adhesive dressing enables firm fixation, minimising device movement or premature displacement



Skin Friendly

Minimises discomfort on application and removal and helps reduce the risk of skin irritation



Waterproof

Outer layer protects site from fluid penetration, allowing patient to shower with the dressing in-situ



ISB403HT

Suitable for:

- Protection and securement of PICC lines
- Protection and securement of Ports
- Protection and securement of PIVC lines
- Protection and securement of CVC lines



ISB202HT

Component details

IV Film Dressing HT Range



ASB301HT



ASB302HT



ASB303HT



ASB304HT




ASB312HT



ASB402HT



ASB403/4HT



ASB202/3HT

REORDER	SIZE	Type	SUITABLE FOR	STERILITY	PCS/UNIT	UNITS/CTN
ASB301HT	5cm x 5.7cm	Bordered	Paediatric, PIVC	Sterile	100	10
ASB302HT	6.5cm x 7cm	Bordered	PIVC	Sterile	100	20
ASB303HT	7cm x 8.5cm	Bordered	PIVC	Sterile	100	18
ASB304HT	7cm x 9cm	Film Only	PIVC	Sterile	100	12
ASB312HT	8.5cm x 10.5cm	Bordered	PIVC, CVC	Sterile	100	16
ASB402HT	6cm x 7cm	Bordered	Paediatric, PIVC	Sterile	100	24
ASB403HT	10cm x 12cm	Bordered	PICC, PORTS, PIVC, CVC	Sterile	100	12
ASB404HT	11cm x 15cm	Bordered	PICC, PORTS, PIVC, CVC	Sterile	100	12
ASB202HT	6cm x 7cm	Film Only	PICC, PORTS, PIVC, CVC	Sterile	50	10
ASB203HT	10cm x 12cm	Film Only	PICC, PORTS, PIVC, CVC	Sterile	25	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.



Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

